

# TiNoI-1<sup>TM</sup>

## NITINOL SPRINGS

PREMIUM NITINOL

by Modern Orthodontics



*Open or Closed*  
*Constant Force*  
*Diamond Polished*  
*Medium Force*





## TiNol-1<sup>TM</sup> Nitinol Springs

**TiNol-1<sup>TM</sup> NiTi open coil spring:** Open coil springs **open and maintain spaces with constant, gentle force**. From start to finish this spring delivers constant, unvarying force for predictable results-even when spring compression is changed the force stays the same.

.010X.030 spring in 15" lengths

**TiNol-1<sup>TM</sup> NiTi open coil spring**

.010" x .030" 15" on spool NT1030O-15I

.012" x .030" 15" on spool NT1230O-15I

**TiNol-1<sup>TM</sup> NiTi closed coil spring with eyelets:**

Use closed coil eyelet springs to close spaces with ease, **without patient discomfort** and within a predictable time period.

.010X.030 spring in 9mm and 12mm lengths in packs of 10

**TiNol-1<sup>TM</sup> NiTi closed coil spring**

.010" x .030" 9mm pk of 10 NT1030C-09M

.010" x .030" 12mm pk of 10 NT1030C-12M

**TiNol-1<sup>TM</sup> NiTi closed coil spring with eyelets for class II correction:**

**Use it in place of class II elastics** when your patient compliance is a problem. Tie it once and it would work permanently for months. Spring should be covered with soft sleeve to eliminate pinching and food trapping.

**TiNol-1<sup>TM</sup> .010X.030 spring in 15 mm and 18 mm lengths in packs of 10**

**TiNol-1<sup>TM</sup> NiTi closed coil spring**

.010" x .030" 15mm pk of 10 NT1030C-15M

.010" x .030" 18mm pk of 10 NT1030C-18M

