

# TiNoI-3F<sup>TM</sup>

## THREE FORCE

PREMIUM NITINOL

by Modern Orthodontics



3 Distinct Force Levels  
Most Physiological Forces  
Diamond Polished  
Laser Marked Midline





## TiNoI-3F™ Progressive Force Nitinol

TiNoI-3F™ Progressive Force wire offers **three force levels** in one amazing archwire

### Anterior Region:



The incisor segment has the **lowest heat-activated force** level specific to moving small, single rooted teeth. The TiNoI-3F™ provides gentle forces for patient comfort, reduces unnecessary tooth pressure and decreases the possibility of root resorption.

### Bicuspid Region:

This segment of the TiNoI-3F™ archwire delivers a **slightly increased force** for moving the larger bicuspids or premolars. The progressively increased force provides quick leveling and alignment to speed up treatment time, without compromising patient comfort levels.

### Molar Region:

The posterior segment of the TiNoI-3F™ generates the **most force** in order to align the larger, more difficult to move first and second molars. Regular archwires have only one force level throughout, which would require 2 or 3 wire changes to achieve the same movement you can get with 1



### TiNoI-3F™ Progressive Force (3Force) Damon Form

● Round	▲ Upper	▼ Lower
.016"	NT16U-DFE3	NT16L-DFE3
.018"	NT18U-DFE3	NT18L-DFE3
● Rectangle	▲ Upper	▼ Lower
.016" X .016"	NT1616U-DFE3	NT1616L-DFE3
.016" X .022"	NT1622U-DFE3	NT1622L-DFE3
.017" X .025"	NT1725U-DFE3	NT1725L-DFE3
.018" X .018"	NT1818U-DFE3	NT1818L-DFE3
.018" X .025"	NT1825U-DFE3	NT1825L-DFE3
.019" X .025"	NT1925U-DFE3	NT1925L-DFE3
.021" X .025"	NT2125U-DFE3	NT2125L-DFE3

### TiNoI-3F™ Progressive Force (3Force) Natural Form

● Round	▲ Upper	▼ Lower
.016"	NT16U-NFE3	NT16L-NFE3
.018"	NT18U-NFE3	NT18L-NFE3
● Rectangle	▲ Upper	▼ Lower
.016" X .016"	NT1616U-NFE3	NT1616L-NFE3
.016" X .022"	NT1622U-NFE3	NT1622L-NFE3
.017" X .025"	NT1725U-NFE3	NT1725L-NFE3
.018" X .018"	NT1818U-NFE3	NT1818L-NFE3
.018" X .025"	NT1825U-NFE3	NT1825L-NFE3
.019" X .025"	NT1925U-NFE3	NT1925L-NFE3
.021" X .025"	NT2125U-NFE3	NT2125L-NFE3

### TiNoI-3F™ Progressive Force (3Force) Euro Form

● Round	▲ Upper	▼ Lower
.016"	NT16U-EFE3	NT16L-EFE3
.018"	NT18U-EFE3	NT18L-EFE3
● Rectangle	▲ Upper	▼ Lower
.016" X .016"	NT1616U-EFE3	NT1616L-EFE3
.016" X .022"	NT1622U-EFE3	NT1622L-EFE3
.017" X .025"	NT1725U-EFE3	NT1725L-EFE3
.018" X .018"	NT1818U-EFE3	NT1818L-EFE3
.018" X .025"	NT1825U-EFE3	NT1825L-EFE3
.019" X .025"	NT1925U-EFE3	NT1925L-EFE3
.021" X .025"	NT2125U-EFE3	NT2125L-EFE3

© Modern Orthodontics Doc#MO/CTG/F01 Issue #1 Rev.00 Dated 25/07/16